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- Behrend EN, Kooistra HS, Nelson R, Reusch CE, Scott-Moncrieff JC. Diagnosis of spontaneous canine hyperadrenocorticism: 2012 ACVIM consensus statement (small animal). *J Vet Intern Med.* 2013;27(6):1292-1304. doi:10.1111/jvim.12192
- Bugbee A, Rucinsky R, Cazabon S, et al. 2023 AAHA Selected Endocrinopathies of Dogs and Cats Guidelines. *J Am Anim Hosp Assoc.* 2023;59(3):113-135. doi:10.5326/JAAHA-MS-7368
- Galac S. Hyperadrenocorticism (Cushing's syndrome) in dogs. In: Ettinger SJ, Feldman EC, Côté E, eds. *Ettinger's Textbook of Veterinary Internal Medicine Expert Consult.* Vol 2. 9th ed. Elsevier; 2024:2004-2021.
- Hess RS. Hypoadrenocorticism. In: Ettinger SJ, Feldman EC, Côté E, eds. *Ettinger's Textbook of Veterinary Internal Medicine Expert Consult.* Vol 2. 9th ed. Elsevier; 2024:2036-2045.
- European Society of Veterinary Endocrinology. Project ALIVE. Accessed June 29, 2025. www.esve.org/alive/intro.aspx
- Bovens C, Tennant K, Reeve J, Murphy KF. Basal serum cortisol concentration as a screening test for hypoadrenocorticism in dogs. *J Vet Intern Med.* 2014;28(5):1541-1545. doi:10.1111/jvim.12415
- Gallego AF, Gow AG, Boag AM. Evaluation of resting cortisol concentration testing in dogs with chronic gastrointestinal signs. *J Vet Intern Med.* 2022;36(2):525-531. doi:10.1111/jvim.16365
- Gold AJ, Langlois DK, Refsal KR. Evaluation of basal serum or plasma cortisol concentrations for the diagnosis of hypoadrenocorticism in dogs. *J Vet Intern Med.* 2016;30(6):1798-1805. doi:10.1111/jvim.14589
- Lennon EM, Boyle TE, Hutchins RG, et al. Use of basal serum or plasma cortisol concentrations to rule out a diagnosis of hypoadrenocorticism in dogs: 123 cases (2000-2005). *JAVMA.* 2007;231(3):413-416. doi:10.2460/javma.231.3.413
- CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples.* 3rd ed. CLSI document EP09c. Clinical and Laboratory Standards Institute; 2018.
- CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline.* 3rd ed. CLSI document EP05 A3. Clinical and Laboratory Standards Institute; October 2014; reaffirmed September 2019.
- CLSI. *Interference Testing in Clinical Chemistry.* 3rd ed. CLSI document EP07 Ed3. Clinical and Laboratory Standards Institute; April 30, 2018; reaffirmed October 2022.

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Catalyst, Catalyst Dx, and Catalyst One are trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

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© 2025 IDEXX Laboratories, Inc. 無断複写・転載を禁止します。 • 09-2692099-00
All other product and company names and logos are trademarks or registered trademarks of their respective holders. IDEXXのプライバシーポリシーは、idexx.comからご確認ください。