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
Images	
	
Cytology Source:	LYMPH_NODE
Clinical History:	Generalized Lymphadenopathy
Pathologist's Report	INTERPRETATION: Limited sample; inconclusive sample- see comments

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📄 🔄 Glucose	97	63 - 114 mg/dL	
📄 🔄 IDEXX SDMA	b 25	0 - 14 μg/dL	
📄 🔄 Creatinine	1.0	0.5 - 1.5 mg/dL	
📄 🔄 BUN	13	9 - 31 mg/dL	
📄 🔄 BUN: Creatinine ratio	13.0		
📄 🔄 Phosphorus	3.3	2.5 - 6.1 mg/dL	
📄 🔄 Calcium	10.0	8.4 - 11.8 mg/dL	
📄 🔄 Sodium	148	142 - 152 mmol/L	
📄 🔄 Potassium	4.5	4.0 - 5.4 mmol/L	
📄 🔄 Na: K Ratio	33	28 - 37	
📄 🔄 Chloride	115	108 - 119 mmol/L	
📄 🔄 TCO2 (Bicarbonate)	21	13 - 27 mmol/L	
📄 🔄 Anion Gap	17	11 - 26 mmol/L	
📄 🔄 Total Protein	5.7	5.5 - 7.5 g/dL	
📄 🔄 Albumin	2.8	2.7 - 3.9 g/dL	
📄 🔄 Globulin	2.9	2.4 - 4.0 g/dL	
📄 🔄 Albumin: Globulin Ratio	1.0	0.7 - 1.5	
📄 🔄 ALT	25	18 - 121 U/L	
📄 🔄 AST	31	16 - 55 U/L	

IDEXX Cancer Dx™

📄 Cancer Dx Lymphoma	a Consistent with lymphoma
📄 Phenotype:	b B cell
a. In dogs with clinical signs concerning for lymphoma, a positive result is supportive of a diagnosis of lymphoma. For dogs under the age of two, confirmation of a positive Cancer Dx lymphoma result should be considered in the absence of strong clinical indications of lymphoma. In dogs without clinical signs being screened for lymphoma due to risk profiles, a positive result suggests an increased likelihood of developing clinical lymphoma. If not performed, a physical exam with lymph node palpation, full history, CBC and comprehensive chemistry testing are indicated. Follow-up in 4-8 weeks to repeat physical exam and lymph node palpation and evaluation for signs of lymphoma is recommended. In rare cases, positive results may indicate other lymphoproliferative cancers, including leukemia and myeloma-related diseases due to shared origin of cells. Results should be interpreted in light of clinical presentation and other laboratory testing in patients with a suspicion of other lymphoproliferative diseases. To learn more about the IDEXX Cancer Dx Test, including recommended next steps and resources, please visit www.idexx.com/CancerDx b. A result consistent with B-cell lymphoma indicates a cancer of B-cell lineage, the most common type of lymphoma diagnosed in dogs. Most B-cell lymphomas carry an improved prognosis compared to T-cell lymphoma.	

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
Images	
Cytology Source:	LYMPH_NODE
Clinical History:	Generalized lymphadenopathy-this is a resubmission (previous cytology was not digital). Tick neg. Normal Blood work
Pathologist's Report	INTERPRETATION: Large cell lymphoma

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1. Vail DM, Pinkerton M, Young KM. Hematopoietic tumors. In: Vail DM, Thamm DH, Liptak JM, eds. *Withrow & MacEwen's Small Animal Clinical Oncology*. 6th ed. Saunders; 2020:688–772. doi:10.1016/B978-0-323-59496-7.00033-5
2. Data on file at IDEXX Laboratories, Inc. Westbrook, Maine USA: Data based on testing performed at IDEXX Reference Laboratories in North America between November 1, 2024, and December 6, 2024. *Analysis Report: IDEXX Cancer Dx Validation, 100282* [008_CancerDx-Validation-Report-2.Rmd].
3. Bailey DB. Hematopoietic tumors. In: Ettinger SJ, Feldman EC, Cote E, eds. *Ettinger's Textbook of Veterinary Internal Medicine*. 9th ed. Elsevier; 2024:2240–2254.